

**designsafe Report**

Application: Bloodless Glucose Monitor  
 Description:

Analyst Name(s): Tom Zhou  
 Company:

Product Identifier:  
 Assessment  
 Type: Detailed  
 Limits:

Facility Location:

Sources:

Risk Scoring System: ANSI B11.0 (TR3) Two Factor

Guide sentence: When doing [task], the [user] could be injured by the [hazard] due to the [failure mode].

Item Id	User / Task	Hazard / Failure Mode	Initial Assessment			Final Assessment		Status / Responsible Risk Level /Reference
			Severity Probability	Risk Level	Risk Reduction Methods	Severity /Comments	Probability	
1-1-1	Doctor conduct tests	electrical / electronic : improper wiring	Moderate Unlikely	Low	Protect circuit from physical contact	Moderate		
1-1-2	Doctor conduct tests	lasers : eye exposure	Serious Unlikely	Medium	Add warning label	Serious		
1-2-1	Doctor calibrate / test / troubleshoot medical monitoring equipment	electrical / electronic : improper wiring	Moderate Unlikely	Low	Protect circuit from physical contact	Moderate		

1-2-2	Doctor	lasers : eye exposure	Serious	Medium	Add warning label	Serious
		calibrate / test / troubleshoot medical monitoring equipment	Unlikely			
2-1-1	Patient	electrical / electronic : improper wiring	Moderate	Low	Protect circuit from physical contact	Moderate
	All Activities		Unlikely			
2-1-2	Patient	lasers : eye exposure	Serious	Medium	Add warning label	Serious
	All Activities		Unlikely			

Privileged and Confidential Information



